


mis<sup>®</sup>  
4BONE<sup>™</sup>  
RCM

Resorbable  
Collagen  
Membrane





## **MIS Bone Grafting Solutions.**


MIS offers a variety of advanced bone and tissue regeneration products aimed at a wide variety of clinical conditions and surgical requirements. All products in this category are packed in varying weights and volumes and in user-friendly containers and utensils, allowing practitioners to choose the precise quantity required for each procedure.

Published by MIS, which reserves the right to ameliorate the products described in this manual as well as to revise this publication at any time and without informing any person of such revision or change. All rights reserved. No part of this publication may be reproduced, transcribed, stored in an electronic retrieval system translated to any language or computer language, or be transmitted in any form whatsoever without the written consent of the publisher.


6.

## Product Description.

4BONE RCM is a resorbable dental membrane made from porcine skin-sourced collagen, aimed for guided bone and tissue regeneration (GBR and GTR). Due to its substantially proven biocompatibility and low immunogenicity, porcine collagen has been successfully used in a variety of medical and dental procedures for many years. 4BONE RCM production involves a cross-linkage technique, enabling functionality as an efficient barrier for a 4-6 months period.



Microporosity 100µm



Microporosity <5µm

### Purified collagen properties

Purified skin-sourced collagen ■ Collagen types I and III ■ Preservation of the fibrous structure (supporting mechanical strength)


### Cross-linkage production technique

A bath in a cross-linking reagent ■ Prolonged resorption time  
■ Neutralization: elimination of the residual cross-linker ■ Chemical analysis of residual reagent

Resorption Control

## Unique Cross-Linkage Technique.

The controlled process of collagen cross-linking allows sufficient time for bone regeneration and rehabilitation of osseous defects. The collagen cross-linking is directly controlled to assure the material's ideal longevity and degradation profile. Verified by animal testing, the degradation of 4BONE RCM membranes is perfectly consistent.



Cross-linkage




Collagen structure

8.

## In-Vitro Resorption Control.

The following comparison demonstrates the rate of mass loss after 48 hours of in-vitro degradation in different types of membranes.

Averages loss of mass



## Advantages.

The 4BONE RCM membrane is the result of extensive research and development, aiming at a functional barrier that remains intact for a specific period of time.

### Efficient

Allowing sufficient time for osseous defects to achieve optimal regeneration, 4BONE RCM barrier enables successful application for various indications.



### Easy handling

Easy to handle, cut to size, shape and apply. Easy to place, with no need for sutures or pins. Flexible and adaptable to various bone topographies.



### Safe

The porcine source ensures an improved sanitary safety and a complete and validated traceability.



### Architectural bone remodeling


The 4BONE RCM microstructure and macrostructure represents a dynamic process, including physico-chemical processes, crystal/protein interactions, cell and tissue colonization and bone remodeling.



10.

#### **4BONE RCM Indications.**

4BONE RCM is a bioabsorbable membrane. It is intended for use in periodontal/dental surgery procedures, for placement in areas of periodontal defect, dental implant, bone defect or ridge reconstruction; to support wound healing post surgery. It is recommended for use with a bone graft material to promote new bone healing.



Tissue preservation is essential for efficient primary closure of the wound and for the correct positioning of the flaps.



## Packaging.

4BONE RCM packaging is designed to ensure sterility. The external layer of the double pouch should be carefully opened and the inner pouch placed onto a sterile field. The membrane should be removed from the inner pouch using sterile gloves or instruments.



4BONE RCM is available in three sizes:


BR-C1525 - 15x25x0.3mm

BR-C2030 - 20x30x0.3mm

BR-C3040 - 30x40x0.3mm









[www.mis-implants.com](http://www.mis-implants.com)

MIS Quality System complies with international quality standards: ISO 13485:2003 - Quality Management System for Medical Devices, ISO 9001:2008 - Quality Management System and CE Directive for Medical Devices 93/42/EEC. MIS products are cleared for marketing in the USA and CE approved.