MIS V3 | 2017




The V3 implant system is the outcome of an exceptionally high-level R&D process that has resulted in an implant that is simple, easy-to-use and offers enhanced functionality and performance. The V3 conical connection implant features built-in design characteristics that provide biological benefits for hard and soft tissue and promotes esthetic results.

\$172C@











© MIS Implants Technologies Ltd. All Rights Reserved.

MIS reserves the right to modify the products described in this manual as well as to revise this publication at any time and without informing any person of such revision or change. All rights reserved. No part of this publication may be reproduced, transcribed, stored in an electronic retrieval system, translated to any language or computer language, or transmitted in any form whatsoever without the written consent of the publisher. Questions, comments or requests will be addressed promptly by contacting MIS specialists directly through our e-mailing address: service@mis-implants.com. This online site highlights current products and reflects all new discoveries and developments.

Note: This user manual is for educational purposes only.

MIS Quality System complies with international

quality standards: ISO 13485:2003 - Quality

Management System for Medical Devices,

ISO 9001: 2015 - Quality Management

System and Medical Device Directive

93/42/EEC. MIS products are CE marked.



Overview

- 8. Introduction
- 9. Raw Material
- 12. Manufacturing Process
- 13. Implant Surface
- 16. Histology
- 17. Hydrophilicity



MIS is a dynamic, state-of-the-art production company, developing and manufacturing a comprehensive range of dental implants designed to provide long-lasting successful solutions to partial and complete edentulous conditions. MIS implant systems combine several advantageous elements such as choice of raw materials, macro-structure, micro-structure and surface treatments, in order to achieve high primary stability and successful osseointegration.

MIS upholds high quality standards by conducting comprehensive quality assurance evaluations throughout the entire production process. The unique MIS implant surface treatment combines sand-blasting and acid-etching to increase surface area, creating both micro and nano-structures and eliminating surface contaminants. The implant surface is continuously monitored by a comprehensive series of tests, conducted both in-house and by internationally recognized research institutes.

Tests include:


- Mechanical tests
- XPS analysis
- Roughness analysis
- Surface analysis
- SEM evauations
- Cytotoxicity tests
- Sterility tests
- Removal torque values
- Histology
- Packaging integrity test




- Biocompatible
- Safe
- Long-term proven clinical success
- Superior mechanical properties

All MIS implants are made from Ti-6AI-4V ELI (Grade 23), the higher purity version of Ti-6AI-4V. This specific type of alloy combines biocompatibility, excellent fatigue strength and low elastic modulus. These benefits make Ti-6AI-4V ELI mechanically superior to titanium grade 4 and the ultimate dental and medical titanium grade.


Similar to commercially pure titanium (Grades 1-4), the outer surface of all MIS implants is comprised of a thin layer of pure titanium dioxide (TiO₂). In this way, bone cells cannot differentiate between the different titanium grades. The TiO₂ layer also prevents metallic ions leaking from the alloy, for safe long-term use.























Roughness (Micro and Nano Structures)


The combination of sand-blasting and acid-etching induces micro and nano-structures that significantly increase surface area of the implant body for optimal osseointegration. The roughened surface improves bone adhesion, as well as the proliferation and differentiation of osteoblasts.



Osseointegration is defined as the attachment of bone to dental implants, and is the critical factor related to the long-term success of dental implants. Osseointegration is determined by both the raw material of the implant, morphology and surface chemical composition.







SEM image of the implant surface

Macro-structure

The geometric design of the body and thread profile of the implant act to increase primary stability and to distribute forces from the implant to the surrounding bone.

Micro and nano-structure

All MIS implants are sand-blasted and acidetched. This surface treatment increases the implant surface area, creating both micro and nano-structures, while eliminating various surface contaminants.

MIS is one of only a handful of companies worldwide using electron microscopy on a daily basis for implant quality inspection.


Sand-blasted and acid-etched surfaces have been substantially proven to maximize the BIC (Bone-to-Implant Contact), achieving rapid and long-lasting osseointegration.

Surface composition


The outer surface of MIS implants, consist of a thin layer of pure titanium dioxide (TiO2). Acid-etching and packaging processes are performed in a controlled environment clean-room to ensure purity and quality. Implants are inspected by scanning electron microscope (SEM) and X-ray photoelectron spectroscopy (XPS), to ensure implants are free of contaminants.



SEM image of the implant surface showing the micro-structure




SEM image of the implant surface showing the nano-structure






Histologic section of a V3 implant, 8 weeks after placement. Courtesy of Prof. Rompen & Prof. Lambert, University of Liege, Belgium.





Current literature demonstrates a linkage between improved bone healing and early osseointegration with the hydrophilicity of surface. MIS implant surface treatment combines sand-blasting and acid-etching. MIS surface treatment ensures surface purity and hydrophilic properties. The images below, demonstrate liquid "climbing" upwards on the implant surface.





Implants


- 18. Introduction: V3 Implant
- 19. Fixture Technical Info.
- 20. External Design
- 21. Implant Range
- 22. Conical Connection
- 24. Procedure



MIS is proud to introduce the V3, an addition to our implant selection, and part of the novel VCONCEPT. The V3 implant features a unique combination of attributes, which result in an innovative implant that provides high initial stability and a state of the art conical connection which incorporates platform switching technology. A large variety of superstructures and components are available, providing solutions for every possible clinical scenario. All superstructures, as are implants, are color coded, according to their restorative platform, with a golden anodized hue for best esthetic results.







- 7 Surface treatment
- 8 Micro-rings



Platform switching

The V3 features platform switching, which keeps the implant-abutment connection away from the bone; minimizing bone resorption. Additionally, platform switching allows more vital growth of the soft tissue.

Conical shape

- The conical root shape of the V3 implant and its unique thread design ensure superior primary stability, making it the implant of choice for a wide range of clinical cases and loading protocols.
- The root shape design makes the V3 an ideal implant when space is restricted due to adjacent teeth or implants.

Three spiral channels and flat apex

- A flat cutting apex allows for final adjustments during placement procedures.
- Three cutting blades at the implant apex establish the self-tapping properties of the V3; supporting a simpler, safer and faster procedure.

Dual thread

 The V3 features a dual thread design which increases the BIC (Bone to Implant Contact) over the entire body of the implant. The dual thread influences implant insertion rate (1.60mm), facilitating a more controlled and faster implant placement.

- The thread profile is especially designed for flawless, easy insertion and high primary stability.
- The V3 is self-tapping with mild bone compression, which enhances primary stability.

Surface treatments

V3 implants are sand-blasted and acid-etched. These surface treatments increase the implant surface area by creating both micro and nanostructures and eliminating various surface contaminants. These treatments ensure surface purity and hydrophilic properties.

Micro-rings

At the neck of the V3, "micro-rings" significantly increase the BIC (Bone to Implant Contact), avoiding bone resorption at the crestal zone.




Length	8mm	10mm	11.50mm	13mm	16mm
Type Ø 3.30mm Screw type implant Narrow platform		V3-10330	V3-11330	V3-13330	V3-16330
Type Ø 3.90mm Screw type implant Standard platform	V3-08390	V3-10390	V3-11390	V3-13390	V3-16390
Type Ø 4.30mm Screw type implant Standard platform	V3-08430	V3-10430	V3-11430	V3-13430	V3-16430
Type Ø 5mm Screw type implant Standard platform	V3-08500	V3-10500	V3-11500	V3-13500	V3-16500

* Implant package also includes: cover screw and final drill.



The V3 features a 12-degree conical connection to ensure a secure fit between the abutment and implant. By minimizing movement at this junction, bone loss is reduced at the crestal level. There is a three-position cone index within the conical connection to help orient the implant during insertion. The cone index also allows for proper abutment positioning. The narrow platform implants include a 3 slot index, while the standard platform implants include a 6 slot index.

STUTIE V






Recommended insertion torque: 35-60 Ncm.








Recommended insertion torque: 35-60 Ncm.







Surgical Procedures

For MIS Implants

- 30. Indications & Contraindications
- 32. Step-by-Step Protocol

Surgical Procedures

Indications & Contraindications

Indications

MIS V3 conical connection implants are intended to be surgically placed in the bone of the upper or lower jaw arches to provide support for prosthetic devices, such as artificial teeth, in order to restore masticatory function. Using a one-stage surgical procedure, the implant allows immediate implantation and immediate function, when good primary stability is achieved and the occlusal load is appropriate.

Narrow implants (Ø3.3mm) are indicated for use in surgical and restorative applications for placement only in the mandibular central, lateral incisor and maxillary lateral incisor regions of partially edentulous jaws, to provide support for prosthetic devices such as artificial teeth, in order to restore masticatory function. Mandibular central and lateral incisors must be splinted if using two or more narrow implants adjacent to one another.



Contraindications

The contraindications customary in oral surgery with other implant materials should

be observed. These include patients taking corticosteroids or anticonvulsants and those receiving radiation or other immunosuppressive therapies. Patients with abnormal laboratory values for BUN, creatinine or serum calcium, patients with diabetes, cardiovascular disease, hypertension above 170/110mm Hq., osteoporotic crush fractures, respiratory disease, thyroid or parathyroid disease, should be excluded as well as patients with diagnosed malignancy during the past five years and those with nodular enlargements, tenderness or an unexplained lump in the head or neck. Implant procedures should not be performed on patients with active osteolythic, or an inflammatory or infectious process in the implant site.



Other Contraindications

- Debilitating or uncontrolled disease
- Hemophilia, Granulocytopenia or other bleeding problems, Steroid use, Prophylactic antibiotics, Brittle diabetes
- Ehler-Danlos syndrome

 Osteoradionecrosis, Renal failure, organ transplantation, Anticoagulation therapy, Idiopathic hypersensitivity, Fibrous dysplasia

Regional enteritis:

- Lack of adequate training of practitioner
- Conditions, diseases, or treatment that severely compromise healing, e.g. including radiation therapy
- Poor patient motivation
- Psychiatric disorders that interfere with the patient's understanding and compliance with necessary procedures
- Unrealistic patient expectations
- Unattainable prosthodontic reconstruction
- Inability of patient to manage oral hygiene
- Patient hypersensitivity to specific components of the implants



Risks

Risks associated with surgical procedures fall into four broad categories:

- 1. Immediate anesthetic and surgical risks.
- 2. Psychological and psychiatric risks.
- 3. Medical threats to long-term retention.
- 4. Long-term deleterious health effects of implants..

The risks may include:

Inadvertent perforation of the nasal and maxillary sinus, local and systemic infections, perforation of soft tissue spaces and nerve injury. Temporary conditions that may result from implant placement may include pain and swelling, speech problems and gingivitis. Long-term problems may include nerve damage, local or systemic bacterial infections, and infectious endocarditis in susceptible individuals. Existing natural dentition may be compromised by improper implant placement. The following list of organ systems with corresponding pathophysiological problems may influence risks:

a) Cardiovascular:

Coronary artery disease, arrhythmias

- b) Respiratory: Chronic obstructive pulmonary disease
- c) Gastrointestinal: Hepatitis, Malabsorption, Inflammatory bowel disease
- d) Genitourinary: Chronic renal failure
- e) **Endocrine:** Diabetes, Thyroid disease, Pituitary/Adrenal disorders
- f) Hematological: Anemia, Leukemia, Blood clotting disorders
- g) Musculoskeletal: Arthritis, Osteoporosis
- h) Neurologic: Stroke, Palsy, Mental retardation



Important Warning

Practitioner's lack of adaquate training, knowledge and experience are considered major risk factors to the patient's health and to the implant's success. Therefore, no implant placement procedure should be performed without prior training by a certified institution.

32. Surgical Procedures Step-by-Step Protocol

This surgical manual is designed to provide an overview of the pre-surgical and surgical procedures applicable to the V3 implant range. Successful implant placement procedures are the result of a wide range of factors. This step-by-step protocol aims to ensure that significant factors are not overlooked.



Step 1.

Patient Selection and Medical History (General medical history)

Patients must be carefully assessed for their ability to safely undergo surgical procedures. Medical history should be evaluated to ensure that patients are not put at risk. Certain medical conditions are considered either absolute or relative contraindications for surgery. These may relate (but not be limited) to the following conditions: Patients who are either taking or have taken medications for the treatment of osteoporosis; immunodeficiency or immunosuppressive treatments; malignancies; head and neck radiation; poorly controlled diabetes or other hormonal disorders; bleeding disorders or anticoagulant therapy; recent myocardia infarction, severe cardiac insufficiency and valve pathology; general bone diseases; hypersensitivity or known allergy to specific relevant materials; psychiatric or personality disorders that limit or interfere with patient understanding and compliance. Please be aware of the fact that updates based on current medical literature may include or exclude certain conditions.



Step 2.

Dental Conditions and Oral Hygiene

A complete and thorough intraoral examination must be performed and recorded. This must include an evaluation of the dentition, oral hygiene, smoking, habits, attitude to oral health and any other relevant information. Implant procedures should not be performed on patients with active osteolitic conditions, active periodontal disease or infectious areas at the implant site. Extreme bruxing and clenching should be taken into consideration.



Step 3. Radiographs and Imaging

Diagnosis and treatment planning for implant placement require the use of different types of radiographs and imaging technologies. Panoramic radiographs are considered standard pre-surgery radiographs, however additional imaging modalities such as CT (Computerized Tomography), tomography and periapical radiographs may be required.

It should be emphasized that certain countries require specific radiographs to be taken before, during and after surgery. It is the obligation of the surgeon to ensure that all required documentation is available and recorded before and after surgery. Vertical and horizontal dimensions of implant sites should be measured and charted. The anatomical relationships of neighboring teeth and proximity to anatomical structures such as the mandibular canal, maxillary sinus and base of the nose must be evaluated. Bone inclination and shape should also be taken into account. Surgical guides with radiopaque markers are recommended. These, coupled with computerized tomographic radiographs may later be altered to be used as computer based surgical guides.



Step 4. Treatment Plan (Patient cooperation)

Based on patient needs, alternative treatment plans should be considered and discussed. The chosen treatment plan should result in a sequence of actions related to initial preparations, the surgical phase and the restorative phase. Surgical Procedures Step-by-Step

Protocol



Step 5A. Implant Selection

V3 implants feature a range of diameters and lengths. V3 Standard platform implants are used in the premolar and molar areas, as well as in the anterior areas. Specific analysis of available bone and distance from vital structures at each proposed site may lead to the choice of specific implant length and diameter; however, current augmentation procedures may allow the use of longer or wider implants.



Each MIS implant comes with labels including all relevant data related to the implant. It is critical that one label is kept as part of the patient's record for future reference.

Step 5B.

Surgical Phase

Surgery should be performed under strict infection control conditions. Preoperative medications and/or antibiotics may be required based on the patient's condition and the extent of surgery and should be decided upon by the operating surgeon. Other monitoring measures, including blood-pressure and pulse measurements should also be considered. An emergency resuscitation apparatus should be available.



Warnings: V3 implants are supplied in a sealed and sterilized package. Implants should never be reused. Implants whose sterility is compromised should not be used. Implants should not be used later than the specific expiration date printed on the package. Implant placement should be performed in accordance with acceptable placement

and loading protocols. MIS recommended procedures are described on pages 24-25. However, it should be emphasized, that procedures recommended by MIS cannot replace the judgment and professional experience of the surgeon.

Step 6.

Osseointegration Phase

Current literature supports multiple loading options. The dentist should decide when to load implants based on specific parameters, related to their individual case.

Step 7. Restorative Phase

V3 implants can support different types of final restorations. Following the solution specified in the treatment plan, the final restoration is fabricated based on accepted restorative protocols. Special attention should be given to ensure correct occlusal adjustment, in order to prevent overloading the implant. MIS superstructures and components must be used with all MIS implants.



radiographs, are recommended. Special attention should be put on oral hygiene and habits, occlusion adjustments and the stability of the prosthesis.



Surgical Kits

- 38. Surgical Kit Description
- 40. Advanced Surgical Instrument Kit
- 42. Kit Contents

The Surgical Kit Surgical Kit Description

The innovative V3 surgical kit is designed for simple and safe implant placement procedures. The kit presents a novel ergonomic design that follows the surgical drilling sequence. In addition, the kit includes a set of length-based pilot drills and color-coded visual cues of both implant diameter and restorative platforms.

The kit, including a surgical ratchet is available for purchaseitem number - MK-T051



Warning! Avoid damage!

The sterilization kit-box and insert must be cleaned and sterilized before each use.


Please see sterilization instructions on page 60.






38.







C1 conical connection insertion tools are supplied separately, MK-0054.



The Surgical Kit Kit Contents

The V3 Surgical Kit includes tools that are designed especially for the step-by-step implant placement process. Correct preparation of the implant site ensures efficient and accurate installation and high primary stability.















Dimensions Material


MT-SMD10	I	Spade marking drill	Length 30mm	Stainless steel
MT-PDM24	I	Position drill mill Ø2.4mm	Ø2.40mm Length 32mm	Stainless steel
MT-PD440	I	Position drill Ø4mm	Ø4mm Length 32.7mm	Stainless steel
CT-P2408	Ι	Pilot drill with built in stopper Ø2.4/2.0 height 8mm	Ø2.40mm Length 28mm	Stainless steel
CT-P2410	Ι	Pilot drill with built in stopper Ø2.4/2.0 height 10mm	Ø2.40mm Length 28mm	Stainless steel
CT-P2411	Ι	Pilot drill with built in stopper Ø2.4/2.0 height 11.5mm	Ø2.40mm Length 33mm	Stainless steel
CT-P2413	I	Pilot drill with built in stopper Ø2.4/2.0 height 13mm	Ø2.40mm Length 33mm	Stainless steel
CT-P2416	I	Pilot drill Ø2.4/2.0 height 16mm	Ø2.40mm Length 37.5mm	Stainless steel
CT-BTC24	Ι	Body try in 2.40mm for tapered impl. procedure	Ø3.20mm Length 28.5mm	Titanium
CT-TDC30	I	Step drill Ø3/2.4mm, external irrigation (red)	Ø2.4-3mm Length 37.5mm	Stainless steel
CT-BTC30	Ι	Body try in 3mm for tapered impl. procedure	Ø3mm Length 28.5mm	Titanium
CT-TDC35	I	Step drill Ø3.5/3mm, external irrigation (blue)	Ø3.5-3mm Length 37.5mm	Stainless steel
CT-BTC35	I	Body try-in 3.50mm for V3 tapered impl. procedure	Ø3.5mm Length 28.5mm	Titanium
CT-TDC40	I	Step drill 4/3mm, external irrigation	Ø4-3mm Length 37.5mm	Stainless steel
CT-BTC40	I	Body try-in 4mm for V3 tapered impl. procedure	Ø4mm Length 28.5mm	Titanium






 CT-NLM30	I	V3 coni. con. long insertion tool for motor, NP	Length 29mm	Stainless steel
 CT-NSM30	I	V3 coni. con. short insertion tool for motor, NP	Length 25mm	Stainless steel
CT-NLR30	I	V3 coni. con. long insertion tool for ratchet, NP	Length 32.4mm	Stainless steel
CT-NLR30	I	V3 coni. con. short insertion tool for ratchet, NP	Length 25mm	Stainless steel
CT-SLM30	I	V3 coni. con. long insertion tool for motor, SP	Length 30.29mm	Stainless steel
CT-SSM30	I	V3 coni. con. short insertion tool for motor, SP	Length 25mm	Stainless steel
CT-SLR30	I	V3 coni. con. long insertion tool for ratchet, SP	Length 32.4mm	Stainless steel
CT-SSR30	I	V3 coni. con. short insertion tool for ratchet, SP	Length 22.3mm	Stainless steel

Dimensions Material





- 50. Using MIS Drills
- 52. Color Code
- 54. Drilling Depth
- 56. Drill Indications
- 58. Final Drill
- 60. Drill Maintenance



Implant placement procedures require the use of several drills with different diameters and characteristics. MIS offers drills with internal and external irrigation, as well as conical and ceramic drills. Most MIS drills are marked for depth control and are color-coded for immediate identification of drill diameter.

Features

MIS drills are available with or sharp drills will ensure atraumatic without internal irrigation. Short drilling procedures and minimal drills are also available for each heat generation. diameter All drills are color-coded The drills are marked for depths of 6, 8, 10, 11.5, 13 and 16mm, and are equipped with a ledge that allows the connection of MIS drill stoppers. All MIS drills have a 120°C cutting degree. The sharpness and high quality of the drills allow for up to 30 uses. Careful use of 6mm 8mm 10mm 11,5mm 13mm 16mm

Drill Stoppers

MIS offers drill stoppers to enable simple and accurate depth control. The C1 Drill Stopper Kits (MK-CDS08, MK-CDS10, MK-CDS11, MK-CDS13) are a series of kits, each used for one specific implant length: 8, 10, 11.5 or 13mm.

For commonly used 3.75 or 4.2 implants, MIS offers a single assorted kit - the C1 Drill Stoppers Kit Standard Platform (MK-BC101), which includes all stoppers required for safe placement of Standard platform implants. All C1 Drill Stoppers Kits are compatible for use with the V3 implant.

C1 Drill Stoppers Kit

C1 Drill Stoppers Kit Standard Platform (MK-BC101)



Diameter	Length
Ø3mm	37.5mm
Ø3.50mm	37.5mm
Ø4mm	37.5mm
Ø4.50mm	37.5mm





Color-coding is used for easy identification of drill or implant diameters as follows:






Standard




Standard

Green

Wide






Depth verification

Depth verification be done by using the Body Try-In tools (CT-BTCxx). (pic.1)

A unique way to estimate the successive implant diameter and the required inter-implant space.

Prior to insertion of a dental implant– the evaluation of the successive implant diameter and the required biological space is a necessity. When coming to evaluate these two parameters, the CT-BTCxx system suggests a unique method – even when only pilot drills have been used and a required correction of the drill location may still be amended. The new suggested method may be used when inserting an implant is required next to a single tooth, between 2 teeth or next to another osteotomy, indicating 1.5mm on each side.

The compatible successive implant diameter is also indicated, as shown in the illustration below. (pic 2)




inimal inter-implant distance estimation

Indication of the successive implants diameter.



(pic.2)





Length & Diameter

The position drill mill is features conical blade geometry of up to 2.4mm and a sharp tip. The drill is 32mm in length and its effective length is 9mm. The drill is made of stainless steel.

The position drill has a diameter of Ø4mm and a sharp tip. The position drill is 32.7mm in length and is made of stainless steel.

V3 Pilot Drills come in five different lengths: 8, 10, 11.5, 13 and 16mm. The first four are equipped with a stopper to simplify the drilling procedure.

Step Drills come in a variety of diameters and lengths.

Aim of Use

The position drill mill is used to mark a reference point for subsequent drills. It is especially useful in immediate placement procedures.


The position drill allows visualization of the actual position of the implant at the end of the drilling procedure. The short, sharp drilling head secures the drill on the bone while the 4mm ring seated above the drill head provides an indication as to the final position of the implant.

Pilot Drills are the first invasive drills used for the preparation of a fixture site. The Pilot Drills are length specific to ensure precise drilling depth.

Step Drills are used to widen the osteotomy. They are NOT length specific, and have laser markings for 6, 8, 10, 115, 13 and 16mm implants. The use of stoppers is highly recommended when using Step Drills.



Final Drill for implant diameters



V3 Implant and Drill Measurements

Each V3 implant package contains a sterile, single-use Final Drill. The drills are recommended for use in bone types 1, 2 & 3. Each Final Drill has a predetermined length and diameter, matching the relevant implant shape and dimension ensuring maximum initial stability while preventing pressure on the implant neck. The length-specific final drills also promote a short and safe drilling procedure. The recommended drilling speed is 200-400 Rpm.



Implant and drill measurements



Correct and careful maintenance of MIS drills is extremely important. Damage to drill tips may cause significant impairment of drill function. The following are detailed instructions for proper maintenance.



Cleaning and Sterilization Instructions

Attention: For your own safety, please wear personal protective equipment (gloves, glasses, mask).

Pre-Cleaning:

1. Soak the drills immediately after use in a detergent and disinfecting solution, preferably an enzymatic cleaning solution, with a pH level between 6-9, prepared with lukewarm water for 5 minutes.

2. Scrub the drills under running water with a soft nylon brush to remove any remaining blood or debris.

3. Rinse under tap water (at least 1 min).

4. Place the drills in a kit, support or rack to avoid any contact between instruments.

Cleaning Procedure for Manual Cleaning:

5. Prepare an ultrasonic bath with a cleaning solution at a concentration and temperature specified in the detergent manufacturer's instructions.

6. Immerse the drills completely and activate the bath for at least the recommended time in the detergent manufacturer's instructions.

7. Rinse under tap water (at least 1 min).

Alternative - Automated Cleaning:

8. Place the rack in a washer-disinfector and apply a cleaning procedure according to the manufacturer's recommendations.

9. Dry on a single-use, non-woven cloth or through a drying cycle of washerdisinfector or filtered compressed air. **10.** Inspect the drills and discard those with defects. Repeat cleaning if necessary..

11. Place the drills in a kit and pack in a sterilization pouch.

12. Steam sterilize according to the chart below. Do not exceed the recommended temperature specified.

13. Keep the sterilization packaging in a dry and clean environment.

Cycle type	Pre-vacuum	Gravity displacement
Temperature	132° C / 270° F	135°C/275°F
Exposure	4 min.	10 min.
Drying time	20 min.	30 min.



Recommendations

- Cutting tools should be used for a maximum of 30 uses.

- Distilled water should be used in order to prevent water spots.

- For all metal instruments, it is recommended to use anticorrosive disinfecting and cleaning agents. They should be Aldehyde free and Ethanolamine free.

- Use only autoclaves that meet the requirements of EN 13060 and EN 285.

- Use a validated sterilization procedure according to ISO 17665.


- For automated cleaning procedures, use a washer-disinfector approved according to EN ISO 15883.





Surgical & Prosthetic Tools

- 64. Ratchet Range
- 66. Specialized Surgical Tools
- 68. Specialized Prosthetic Tools
- 72. Screw Tests
- 73. Maintenance



MIS offers a line of uniquely designed ratchets, to simplify both prosthetic screw tightening and implant insertion, allowing for an accurate and safe performance. To prevent damage to the mechanism, it is critical that the ratchet is used only with keys and adapters that are specifically designed for it.

Three ratchet types, to allow an accurate and safe procedure:





Warnings

MIS recommends the use of a torque controlled driver whenever possible.

The ratchet wrench MT-RI030 may transfer torque levels that do not correlate to the recommendations specified for implant placement or screw fastening.

Excess loads may result in damage to implants, components, screws, and even to the bone-to-implant interface.

Beware that the recommended torque for implant placement is 35-60 Ncm.



Instrument Maintenance


- The device is delivered non-sterile.
- Cleaning and sterilization are required prior to use.



Cleaning and Sterilization

For cleaning and sterilization instructions please refer to page 73.

User Instructions:




Store the ratchet on its own, not attached to any tools.

Clean thoroughly, immediately after use.





V3 Insertion Tools

V3 implants are divided into Narrow platform implants (3.30mm) and Standard platform implants (3.90, 4.30 and 5mm). Long and short insertion tools are available for each of the V3 platforms, for hand-piece connection and for use by ratchet/manually.









Friction Fit MT-RE172/ MT-RE160

The friction fit extractors are designed to separate the friction fit abutments from the implant. The extractors are color-coded; purple for Standard abutments and yellow for Narrow abutments.



Mode of Action

The Extractor Key applies vertical load parallel to the long axis of the implant. Thus it can release a "locked" abutment from an implant.


For Standard Implants






For Narrow Implants







SOS Broken Screw Kit

MT-TF172 / MT-RT001/ MT-HW001/ MT-TF160/ MT-RT002

The SOS Broken Screw Kit was designed to facilitate the removal of a broken screw from within an implant.




SOS Tools









Thread Former MT-TF160

Thread Former MT-TF172

Hand Wrench MT-HW001

Retriever MT-RT001

Retriever MT-RT002

Instructions for use:



1.

A. Connect the retriever to a micromotor.

B. Adjust the micromotor to low speed (15-25 RPM), max. torque and in reverse mode.



۷.

A. Apply mild pressure with the retriever to the top of the broken screw.

B. While maintaining pressure, activate the motor. This action should release the screw. If the screw is still not released, apply intermittent pressure on the screw.



З.

If internal threads are damaged:


A. Use the thread former with care.

B. Be sure to align the thread former parallel to the long axis of the implant.

C. Always start by using a hand wrench. Apply gentle but firm force while turning the thread former in a clockwise direction. Release the pressure at the end of each complete turn by turning it 30° in a reverse direction, and repeat the action as needed.

D. In instances where greater torque is needed, a ratchet may be used.





Instrument Maintenance

Attention:

For your own safety, please wear personal protective equipment (gloves, glasses, mask).

Pre-Cleaning:

- Disassemble the device if required.
- Soak all instruments immediately after use in a detergent and disinfecting solution, preferably an erzymatic deaning solution, with a pH level between 6-9, prepared with lukewarm water for 5 minutes.
- Scrub the instruments under running water with a soft nylon brush to remove any remaining blood or debris.
- Rinse under tap water (at least 1 min).
- Place the instruments in a kit, support or rack to avoid any contact between them during the next cleaning procedure.

Cleaning Procedure Manual Cleaning:

- Prepare an ultrasonic bath with a cleaning solution at a concentration and temperature specified in the detergent manufacturer's instructions.
- Immerse the instruments completely and activate the bath for at least the recommended time in the detergent manufacturer's instructions.
- Rinse under tap water (at least 1 min).

Alternative: Automated Cleaning

Place the rack in a washer-disinfector and apply a cleaning procedure according to the manufacturer's instructions.

Drying and Sterilization

- Dry on a single-use non-woven cloth or through a drying cycle of a washer-disinfector or with filtered compressed air.
- Inspect the devices and discard those with defects.

Repeat cleaning if required.

- Assemble the device if required.
- Place the devices in a kit, and pack in sterilization pouch.
- Steam sterilize according to the table below.
 Do not exceed the recommended temperature specified.
- Keep inside the sterilization packaging in a dry and clean environment.

Cycle type	Pre-vacuum	Gravity displacement
Temperature	132°C/270°F	135°C/275°F
Exposure	4 min.	10 min.
Drying time	30 min.	30 min.



Recommendations

- Cutting tools should be used for a maximum of 30 uses.
- Distilled water should be used in order to prevent water stains.
- For all metal instruments, it is recommended to use anticorrosive disinfecting and cleaning agents. They should be Aldehyde, Ethanolamine, chlorine and acid free.
- Use only autoclaves that meet the requirements of EN 13060, EN 285.
- Use a validated sterilization procedure according to ISO 17665.
- For Automated cleaning procedures, use a washer-disinfector approved according to EN ISO 15883.



Packaging

- 76. Implant Package
- 78. Implant Identification Codes
- 80. Implant Package Handling
- 88. Transparency
- 89. Symbols



The innovative MIS packaging system is designed for simple and easy use. All of our implant boxes feature distinctive colors, large typeface, clear data labels and a pull tab for quick opening. Boxes are a uniform shape and height, specifically designed to fit in clinic cabinets for easy accessibility and compact space-saving storage.







A convenient 4-implant package is available. The drawer-like box is ideal for storage in drawers or cabinets for easy identification of implant type, diameter and length.




Implant removal from the sterile inner tube

Double Container System

To ensure that implants are sterile and to prevent surface contamination, each implant is stored in a titanium sleeve within an internal plastic tube. This tube is held in a larger sealed outer tube, marked with all relevant information. The inner tube is therefore sterile, and may be brought into the sterile surgical field whenever needed. An anti-rotation mechanism inside the titanium sleeve insures a safe implant removal.






Identification markings enable quick identification of implant diameter (top), implant length (center) and implant platform size (bottom).



Identification labels:







The distinctive blue V3 implant boxes are a uniform shape and height, specifically designed to fit in clinic cabinets for easy accessibility and compact space-saving storage.



Fig. 1

The convenient pull tab facilitates easy and quick opening of the box during surgery.



Fig. 2

Open the outer tube by turning the cap counter-clockwise. Drop the sterile inner tube into the sterile field.



Fig. 3

The implant is held by the titanium sleeve. To expose the implant - hold the tube with the titanium sleeve facing up. Rotate and pull to open the upper cap.



82.



Use one of the following options to remove the implant from the inner tube:

Fig. 4 Contra-angle hand piece



Fig. 5 Manual and Ratchet connection



Fig. 6 Implant placement (illustrated manually).



Fig. 8

Remove the cover screw using the CT-SLR30/ CT-SSR30 key





Fig. 9

Attach the cover screw to the implant using the CT-S#R30 key. Tighten the cover screw using the MT-RDL/S30 KEY




Fig. 10


Attach the data label from the implant package onto the dental patient records.

84.

Implant Data Label

Each package contains three data labels, which include all required information pertaining to the implant. The following image illustrates the label:











MIS offers a planning transparency, illustrating the full V3 implant range. It includes two sets of images: one at actual size, and the other at a magnification of 125%, relevant for use with panoramic radiographs that include a similar inherent magnification. In addition, the transparency includes a 1:1 ruler.

By placing the appropriate section of the transparency on a radiograph, a clinician may choose the best fitting implant diameter and length, as part of the planning process.

The transparency available for V3 implants is: Cat No. MP-CONV3







Key to symbols on labels and instruction leaflets:









All rights reserved. No part of this publication may be reproduced transcribed, stored in an electronic retrieval system, translated into any language or computer language, or be transmitted in any form whatsoever, without the prior written consent of the publisher. Warning Only a licensed dentist should use these products. MP-UI0V3 Rev.2

MIS Implants Technologies Ltd. www.mis-implants.com

MIS Quality System complies with international quality standards: ISO 13485:2003 – Quality Management System for Medical Devices, ISO 9001: 2015 – Quality Management System and Medical Device Directive 93/42/EEC. MIS products are CE marked.